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Press the button ③ up, and the dental simulation will move in the direction shown by arrow (upwards). When the

Figure 16

required position is reached, release the button and the dental simulation will immediately stop moving (if the button is pressed down and held, the dental simulation will automatically stop when it reaches the extreme position)

D. Lowering button

Press the button ③ down, and the torso of dental simulation will move in the direction shown by arrow (downward). When the required position is reached, release the button and the dental simulation will immediately stop moving (if the button is pressed down and held, the dental simulation will automatically stop when it reaches the extreme position)

E. Backward-tilting button

Press the button ③ to the left, and the torso of dental simulation will move in the direction shown by arrow (backward). When the required position is reached, release the button and the dental simulation will immediately stop moving (if the button is pressed down and held, the dental simulation will automatically stop when it reaches the extreme position)

F. Forward-tilting button

Press the button ③ to the right, and the backrest of dental simulation will move in the direction shown by arrow (forward). When the required position is reached, release the button and the dental simulation will immediately stop moving (if the button is pressed

down and held, the dental simulation will automatically stop when it reaches the extreme position)

4. Water storage bottle

- A. Clean water storage bottle(600ml), the clean water use for handpiece and syringe comes from here.
- B. Waste water storage bottle(600ml), the waste water from the suction will be stored up here.


Figure 17

5. High/ Low speed handpiece

- A. Connect water, air and power supply. Open the general air switch on the side of the movable simulation treatment machine, and check the pressure gauge after open the plastic door . (Figure 18) The value should be 0.5 ~ 0.6 MPa (factory setting). Adjust the filter relief valve if it is required to maintain the said value. Open the plastic door, pull the handle on the top of the filter relief valve up for about 10mm as shown in (Figure 18), turn the handle clockwise to increase the pressure and anticlockwise to decrease the pressure.
- B. Take the handpiece from the holder, step the pedal switch for operation. Be noted that the pressure indicated on the pressure gauge of the instrument disc is the operating pressure of the handpiece, which should be no more than the rated maximum pressure of the handpiece to protect the handpiece against damage (High speed: 2.0-4.0bar, Low speed: 3.5-4.0bar), see the Figure 19. Adjust the operating pressure of handpiece if it is required by regulating the main control valve under the instrument disc. Turn the handle clockwise to increase the pressure and anticlockwise to decrease the pressure. Adjust carefully and slowly.

Figure 18


Figure 19

6. Three-way syringe

The left button is for water and the right one is for gas.

See figure 20.


Figure 20

7. Weak suction

Saliva aspirator is provided with this equipment. Take the saliva ejector from the holder to start operation immediately. Water connection is required for low aspiration to guarantee minimum operating pressure required.

8. Clean water bottle

Water for handpiece is directly from the water bottle, therefore the bottle shall be supplied with medical distilled water on time, with water filling described as below: Turn off the air switch beside the water bottle firstly. After all the compressed air in the bottle discharged, hold the water bottle securely and turn clockwise to take it off. Then fill water in it, turn it on counter-clockwise, until it fixed on the bottle cup tightly (Air tightness must be regarded). Finally turn on the air switch.

9. Waster water bottle

The waste water comes out from the simulation mouth through the saliva aspirator, will go into the waste water storage bottle. (On the left side of the machine without blue tube in it.) Hold the water bottle with both hands, rotate clockwise to take off the water bottle, pour out the waste water and rotate anticlockwise to tighten the bottle (sealed).

● Maintenance

- (1) After adjust the simulation head, ensure it is locked before it is used.
- (2) Regularly cleanse the water filter.
- (3) Power supply is 220VAC.
- (4) Cut the power supply before repair the wearable component and cleanse, maintain the treatment machine.
- (5) Should close the lamp, when it is not used.
- (6) To ensure the neat and tidy of the treatment machine, cleanse the surface of the machine and simulation with hospital use alcohol regularly is suggested.

● Note

- (1) The power cord should be configured as standard and the ground wire should be firmly connected.
- (2) When replacing electronic components, the power must be turned off.
- (3) Before the maintenance and cleaning of the equipment, the power must be turned off.

● Others

1. Transport and storage environment:

- (1) Ambient temperature: $-40 \sim +70$ °C.
- (2) Relative humidity: 20% ~ 90%, including condensation.
- (3) Pressure: 86 ~ 106 kPa.
- (4) The rain must be prevented during transportation and gently handled to avoid vibration.
- (5) Treatment of waste water and other materials must comply with local environmental protection regulations.
- (6) Packaging units should be stored in places where the relative humidity does not

exceed 80%, where there is no corrosive gas and air circulation.

(7) The maintenance of the equipment must be performed by professional technicians designated by our company. If the user disassembles and repairs the device by himself, the device may be damaged, and if this happens, our maintenance service will no longer be available.